

Bath chairs/ Bathing Systems

Below are just some examples of bath chairs/bathing systems. Please use this a reference to search the websites listed to see other options available.

Standard Bath Chair



www.drivemedical.com

- User must have good static and dynamic sitting balance
- No arm rests, must be able to perform sit to stand transfer independently or with assist they have available (if safe)
- Would still need to step over tub to sit in this chair
- Weight Capacity is 400 pounds
- Seat height 15.5" - 19.5"
- Seat dimensions 19.25" (W) x 11.5" (D)
- Other versions of this chair available, with arms rests and padding (more difficult to get insurance to cover)

Standard Transfer Tub Bench



www.drivemedical.com/us/en

- Allows user to sit on bench on the outside of the tub and scoot over to the inside of the tub
- Grab bar can go on either side of bench
- User should still have good static and dynamic sitting balance
- Consider skin integrity with scooting across bench
- Weight Capacity is 400 pounds
- Seat height 17.5" - 22.5"
- Seat dimensions 26" (W) x 18.5" - 19.5" (D)
- Other versions of this bench available, with padding or combination with commode (more difficult to get insurance to cover)

Bath chairs/ Bathing Systems

Rifton WAVE bath system

www.rifton.com/



- Offers FOUR base options
- Free standing (to be placed in the tub)
- Tub stand (raises height of bath chair 5.5" or 11.5")
- Tub transfer base (allows you to easily transfer your child to the bath chair on the outside of the tub, rotate and slide them into the tub. Less caregiver strain)
- Shower stand (for roll in shower)
- Offers multiple accessories and supports
- Chair is easily adjustable
- Comes in three sizes, based on user's height and weight
- Max user height 74 inches and max user weight 250 pounds

Otter Bath Chair

www.inspiredbydrive.com/



- Comes in four different sizes
- Can be used inside of the tub and with a roll in shower
- Max user height 72 inches, Max user weight 250 pounds
- Has adjustable legs to add 7 inches in height
- Can get a wheeled shower base for roll in showers and an additional tub stand
- Comes with multiple accessories for safe bathing

HTS Hygeine and Toileting System

www.rifton.com/



- Can be used for toileting and bathing (only with a roll in shower!)
- Consider this type of chair if the user performs a bowel management program that might require them to sit over a toilet for extended periods.
- Can be attached to a toilet or on a mobility base with or without a commode
- Multiple supports available for head, trunk and LE positioning
- User can be positioned in the most optimal way for toileting.
- Waterproof foam padding
- Max user height 74 inches and max user weight approximately 250 pounds
- Max user weight to utilize foot board for transfers is 150 pounds