

Self-Propelled Manual Wheelchairs

Below are just *some examples* of manual wheelchairs. Please use this a reference to search the websites listed to see other options available.

Very early on

Bella's Bumbas



bellasbumbas.com

www.facebook.com/bellasbumbas

- “Bumbo” type seat with platform and large wheels and hind wheel
- Provides an option for mobility very early on
- Non-profit organization
- Families ONLY pay for shipping costs
- Based out of Webster, NY

Folding Wheelchair

- More easily transported in a vehicle
- These wheelchairs are able to grow with your child in width, depth and height
- These chairs will have swing away footrests and foot plates
- Folds “hotdog style”, need to remove back rest and seat cushion to fold

Zippie XCape

www.sunrisemedical.com



- Weight of wheelchair is approximately 19 pounds
- Seat width 8” - 18”
- Seat depth 8”-20”
- Seat to floor height is 13” - 18.5”
- Camber* adjustment 0 and 3 degrees
*camber: adjustments made to wheels to increase maneuverability and stability

Zippie 2

www.sunrisemedical.com



- Weight of wheelchair approximately 25 pounds
- Seat width 12” - 16”
- Seat depth 12” - 18”
- Seat to floor height is 15” - 20”
- No camber adjustments

Rigid Frame Wheelchairs

- Non-Folding- However, seat backs will flip down and wheels are removable
- These wheelchairs are able to grow with your child in width, depth and height
- Rigid foot plates- they do not swing out of the way
- These chairs are great for recreational use and ideal for children with good trunk control and upper body strength
- Requires less effort to propel compared to folding wheelchairs

Zippie Zone



www.sunrisemedical.com

- Typical weight of wheelchair is pounds 14 pounds
- Seat width *- 16"
- Seat depth 8" - 16"
- Seat to floor height 14" - 19"
- Can adjust camber* of wheels at various degrees
*camber: adjustment made to wheels to increase maneuverability and stability

Little Wave



www.kimobility.com

- Standard weight of wheelchair is 16 pounds
- Seat width 10" - 19"
- Seat Depth 10" - 16"
- Seat to floor height 11.5" - 20.5"
- Can adjust camber* of wheels at various degrees

TiLite Pilot



permobilus.com

- Light weight carbon fiber frame, 11 pounds without rear wheels
- Seat width 8" - 14"
- Seat Depth 9" - 15"
- Seat to floor height 12" - 16.5"
- Offers different wheel lock options than other chairs as well as foot rest adjustments

Wheelchairs requiring caregiver assistance

Adaptive Strollers

Not always recommended to go through insurance for adaptive strollers if a wheelchair is recommended in the near future. Speak with your therapist or vendor to decide the best option.

Special Tomato Push Chair

www.specialtomato.com

www.adaptivemall.com



- Approximate weight of stroller is 22 pounds
- Maximum weight capacity is 90 pounds
- Target user height is 30"—50"
- Ability to be reclined
- Accessories and additional inserts can be added for additional support
- Do not need to remove wheels to fold it

Zippie Voyage

www.sunrisemedical.com



- Approximate weight of stroller with base is 28 pounds
- User weight capacity is 75 pounds
- Reversible seat with tilt and recline
- Indoor base is available for home use
- Multiple accessories are available to be added on to the stroller

Convaid Trekker Stroller

www.convaid.com



- Weight of frame (without seat) is 27 pounds
- Weight capacity up to 110 pounds
- Maximum seat width is 14"
- Tilt and recline seating
- Reverse seating available
- High-low indoor frame can be used to transition child from properly positioned in their wheelchair to sitting at the dining room table at home
- Many accessories available such as a stroller canopy, oxygen tank holder, harness's, straps and supports.

Tilt in Space Wheelchairs

- Provides added support and positioning options ideal for children with challenges with tone or head/ trunk weakness that requires a break from gravity (Tilt)
- Tilt is also beneficial for children requiring assistance with breathing, feeding, performing catheterizations and pressure relief.

Little Wave Flip



www.kimobility.com

- Weight of wheelchair without wheels, foot rests or seating system
- Max weight capacity is 165 pounds
- Will fold after removing seating system
- Seat width 10" - 18"
- Seat depth 12" - 20"
- Multiple accessories available

Zippie Iris



www.sunrisemedical.com

- Weight of wheelchair is approximately 40 pounds
- Has the ability to be folded with removal of seat back and cushion
- Rotation in space technology, so when you tilt the chair back it will not continue to tip over like a standard stroller.
- Seat width 10" - 16"
- Seat depth 10" - 18"